

Instructions for use

Product description


Instrument Disinfectant and Cleaner

Date of revision: 09/03/23 v1.6


topdental™
Bossklein®

Disinfecting and cleaning concentrate for instruments, burs and endodontic instruments and items. Chloride, aldehyde and phenol free with a wide material compatibility and broad microbiological efficacy.


2L and 5L bottle



1. Add 20ml ± 1ml to 980ml of clean water for a 2% solution, 40ml ± 1ml to 960ml for a 4% solution.



2. Completely immerse the instruments into the solution ensuring that no air is left within hollow instruments.




60 mins

3. To achieve full efficacy use a 4% solution for a 60 minute soaking time.




5. After disinfection, rinse instruments thoroughly with RO or distilled water prior to further processing and sterilisation. Follow the reprocessing instructions from the instrument manufacturer. Change the solution daily. For professional use only. See product labels and MSDS for further information.


1L dispensing bottle



1. Using the dosing device on the bottle, add 20ml ± 1ml to 980ml of clean water for a 2% solution, 40ml ± 1ml to 960ml for a 4% solution.



2. Completely immerse the instruments into the solution ensuring that no air is left within hollow instruments.



60 mins

3. To achieve full efficacy use a 4% solution for a 60 minute soaking time.



5. After disinfection, rinse instruments thoroughly with RO or distilled water prior to further processing and sterilisation. Follow the reprocessing instructions from the instrument manufacturer. Change the solution daily. For professional use only. See product labels and MSDS for further information.

Efficacy 2% (15 min): Bactericidal incl. MRSA. Yeastcidal. All enveloped viruses: Vaccinia virus, HBV, HCV, HDV, HIV, Influenza, Coronavirus, Herpesviridae, Filoviridae, Rubella, Rabies, Measles.

Efficacy 4%: As 2% and to include Non-enveloped viruses (30 min): Murine norovirus, Adenovirus 5. Mycobactericidal & Tuberculocidal (60 min).

Tested to Standards: EN 14885, EN 14561, EN 14562, EN 14563, EN 14476, EN 14348, EN 13727, EN 13624.

Composition: Aqua, N-(3-Aminopropyl)-N-Dodecylpropane-1,3-Diamine, N,N-Didecyl-N-Methyl-Poly(Oxyethyl)Ammonium Propionate, Alcohols, C9-11, Ethoxylated, Tetrasodium N,N-Bis(Carboxylatomethyl)-L-Glutamate, Monoethanolamine, Corrosion Inhibitor, Aroma, Dye.

Danger: H315: Causes skin irritation. H318: Causes serious eye damage. H410: Very toxic to aquatic life with long lasting effects. P280: Wear protective gloves/protective clothing/eye protection/face protection. P301+P330+P331: IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. P302+P352: IF ON SKIN: Wash with plenty of water. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **Storage:** P403+P233: Store in a well-ventilated place. Keep container tightly closed. P235: Keep cool. P410: Protect from sunlight.

Contraindications. Conditions to avoid: Heat. **Disposal.** Do not discharge concentrate into drains or rivers.

Intended Users. The products are used within the dental surgery by the dentist or dental assistant as part of their general disinfecting routine. Patients are usually not present during use, but all patients that are seen in the surgery contribute to the need for the products.

Any serious incident related to this medical device must be reported to the manufacturer and the competent authorities.

UFI: 3X00-W0GP-900K-VWW7



1639



DIRECTA AB
P.O. Box 723, Finvids väg 8 SE-194 27
Upplands Väsby, Sweden



For additional languages
see www.topdental.co.uk



Danger



DIRECTA topdental ORSING
parkell KOHLER RÖNVIG
TrollDental Sendoline CONTACTEZ
physics ceramir