

Tetric® PowerFill | Tetric® PowerFlow

1

Self-Etch



OR Etch & Rinse



15-30 s



10-15 s

OR Selective Enamel-Etch



15-30 s



2

Adhese® Universal
VivaPen®



20 s



3

Tetric®
PowerFlow



Optional



3 s

4

Tetric®
PowerFill



3 s



Bluephase®
PowerCure



Direct posterior restorations, light-cured from occlusal aspect
Direkte Seitenzahnrestorationen, okklusale Belichtung

	Tetric PowerFill	Tetric PowerFlow	Adhese Universal
500 – 900 mW/cm ²	20 s	20 s	10 s
900 – 1400 mW/cm ²	10 s	10 s	10 s
1800 – 2200 mW/cm ²	–	5 s	5 s
2700 – 3300 mW/cm ²	3 s	3 s	3 s

see Instructions

Date information prepared:
2021-07-02
693240/Rev. 1

CE 0123

Ivoclar Vivadent AG
FL-9494 Schaan/Liechtenstein

ivoclar
vivadent